TapeTrack Documentation - https://rtfm.tapetrack.com/

# TMSS10LibraryHealth

# Synopsis

TMSS10LibraryHealth [-S logon string] [-V volume filter] [-i ignore no library health option]

# Description

TMSS10Libraryhealth produces a report on each customers media health status.

# Options

- -S The logon string argument sets the Server logon information
- -V Volume Filter which can be used to filter the traversal of Customer, Media, Volume, repository and flag values.
- -i Ignore the **no library health** option.

## Notes

The No Library Health can be set on the media properties options tab.

## **Technical Support**

The TapeTrack Software is commercially supported by a full time help desk staff. If you are experiencing problems or want some advice on how to configure or use the product please see the Accessing Technical Support page.

## **Exit Statuses**

- 1. zero Program has ended successfully.
- 2. non-zero Program has not ended successfully.

## Environment

TMSSSERVERPROXY If defined the program will route all TapeTrack TCP/IP traffic through a HTTPS proxy. The value of the variable should be in the format user:password@host:port. To debug the proxy connection use variable TMSSAPILOGDIR.

TMSSAPILOGDIR If defined the program will write out a trace file to this directory.

TMSSPWPATH When no password value is passed in the logon string the program will look for the password in file C:/tapetrack/pw/batch, where user is the user value passed in the logonstring. If you wish to change this default path, you can set the path in TMSSPWPATH.

TMSSNOMD5 If defined the program will not MD5 hash passwords before sending them to the TapeTrack Server. This is required when relying on Windows Active Directory authentication. It should be noted that although the password is not being hashed, it is still being encrypted during transmission.

#### **Files**

stdout: Output Reports.
stderr: Diagnostic messages.
stdin: Program input.

## Example

<pre>tmss10libraryhealth -S user:-password@server -V "US02.*.*"</pre>		
System: TapeTrack Report Description: Library Health Report		
Report Notes: Customer(Los Angeles Data Center) ì-□□! Report Width: 63 bytes		
Media	Health Status	
LTO Cartridge	Not ideal Reason	Penalty
	No Sync	- 30
	No DR Strategies for 4 days	-10
	No Offsite activity for 4 days	-10
Overdue Audit Warni	Volume-ID Current Target .ng Audit Date	Last-Scanned

From: https://rtfm.tapetrack.com/ - TapeTrack Documentation

Permanent link: https://rtfm.tapetrack.com/cli/tmss10libraryhealth?rev=1506555188

Last update: 2025/01/21 22:07



2025/05/08 14:06