

Example

Command line syntax

Move all US02 LTO volumes with a barcode ending in 5L6.

```
TMSS10MoveLocation -S user:-password@server -V "US02.LT0.*5L6" -f "LIBR" -t  
"OFFS"
```

Output

Seq.	Barcode	From	To
Container			
<hr/>			
1	US02.LT0.000015L6	Library	Offsite Vault
2	US02.LT0.000035L6	Library	Offsite Vault
3	US02.LT0.000065L6	Library	Offsite Vault
4	US02.LT0.000075L6	Library	Offsite Vault
5	US02.LT0.000085L6	Library	Offsite Vault
6	US02.LT0.000095L6	Library	Offsite Vault
7	US02.LT0.000105L6	Library	Offsite Vault
8	US02.LT0.000115L6	Library	Offsite Vault
9	US02.LT0.000125L6	Library	Offsite Vault
10	US02.LT0.000135L6	Library	Offsite Vault
11	US02.LT0.000145L6	Library	Offsite Vault
12	US02.LT0.000155L6	Library	Offsite Vault
13	US02.LT0.000165L6	Library	Offsite Vault
14	US02.LT0.000175L6	Library	Offsite Vault
15	US02.LT0.000185L6	Library	Offsite Vault
16	US02.LT0.000195L6	Library	Offsite Vault
17	US02.LT0.000205L6	Library	Offsite Vault
18	US02.LT0.000215L6	Library	Offsite Vault
19	US02.LT0.000225L6	Library	Offsite Vault
20	US02.LT0.000235L6	Library	Offsite Vault
21	US02.LT0.000245L6	Library	Offsite Vault
22	US02.LT0.000255L6	Library	Offsite Vault
23	US02.LT0.000265L6	Library	Offsite Vault
24	US02.LT0.000275L6	Library	Offsite Vault
25	US02.LT0.000285L6	Library	Offsite Vault

26	US02.LT0.000295L6	Library	Offsite Vault
27	US02.LT0.000305L6	Library	Offsite Vault
28	US02.LT0.000315L6	Library	Offsite Vault

From:

<https://rtfm.tapetrack.com/> - **TapeTrack Documentation**



Permanent link:

<https://rtfm.tapetrack.com/cli/tmss10movelocation?rev=1511747962>

Last update: **2025/01/21 22:07**