Run TMSS10Ping From Lite

Ping TapeTrack Server allows you to ping the TapeTrack Framework Server, using TMSS10Ping, from within Lite, negating the need to open a command prompt and run the program from there.

Unlike the ping command which sends ICMP echo packets, TMSS10Ping sends payload packets via the TCP/IP protocol. This means that even if a firewall is configured to block ping requests, if TapeTrack connections are allowed, then TMSS10Ping packets will be successful.

The settings to execute TMSS10Ping via the Lite menu are contained in the Lite Config file. **These** settings need to be customised to your values before the command will execute successfully. Default settings will connect to localhost.

Execute TMSS10Ping

From the drop down menu in Lite select Ping TapeTrack Server.

Lite will open a command prompt window and execute TMSS10Ping and display the results.

Successful Connection

A successful connection will show four byte packets with the travel time as well as a summary on round trip minimum, average and maximum times.



Click Yes to the Command Complete popup to close the Command Prompt, click No if you want to leave the Command Prompt open.



Unsuccessful Connection

Unknown Server

If the server address has not been resolved to a correct IP address, the error below will be shown.



Contact your Server Administrator to ensure the correct IP address is entered into the Lite config file.

Click Yes to the Command Complete popup to close the Command Prompt, click No if you want to leave the Command Prompt open.



No Server Reply

If the IP address resolves correctly, but no reply is received from the server TMSS10Ping will return a connect to server failed, WSA Error(10061) error.

Run Command		\times
Output		
TMSS10Ping: Connect to TapeTrack Server failed: WSAError(10061)		

Contact your Server Administrator to ensure the TapeTrack Framework Server is running and your firewall is not blocking access to the Server ports.

Click Yes to the Command Complete popup to close the Command Prompt, click No if you want to leave the Command Prompt open.

Commar	nd Complete	×
?	Your command has completed with return code 1. Would you like to close the command output window?	
	Yes No	

From: https://rtfm.tapetrack.com/ - **TapeTrack Documentation**

Permanent link: https://rtfm.tapetrack.com/lite/function_ping?rev=1611626271

Last update: 2025/01/21 22:07



3/3