


# Relaxed Barcode Matching

The [Media Option](#) Relaxed Barcode Matching will, when a scanned Barcode is not found, attempt to match the Barcode with any [Volumes](#) with a L|T|J suffix.



## Options

- True: Enables Relaxed Barcode Matching.
- False: Disables Relaxed Barcode Matching.

# Examples

If a barcode 100123 is scanned and not found, Relaxed Barcode Matching will then attempt to match against the barcodes L|T|J suffix, such as 100123L4, 100234T1 or 100234J2.

[options](#), [media options](#), [barcode](#)

From: <https://rtfm.tapetrack.com/> - **TapeTrack Documentation**

Permanent link: [https://rtfm.tapetrack.com/options/media/relaxed\\_barcode\\_matching?rev=1573428224](https://rtfm.tapetrack.com/options/media/relaxed_barcode_matching?rev=1573428224)

Last update: **2025/01/21 22:07**

