

Labelling Promedia Drawers

Labelling the Promedia Drawers in the following manner allows the configuration of TapeTrack allowing operators to quickly locate a [Slotted Volume](#), or locate a specified [Slot](#) to place a [Volume](#) in.

With TapeTrack [Zones](#) and [Slotting](#) information configured correctly, each [Slot](#) will be referenced with the [Zone-ID](#), Drawer row and [Slot](#) number in the format of [Zone-ID](#) Row#.Slot#.

For example **Red-2 3.16** would be in the 2nd drawer (Red-2), in the third row, 16th [Slot](#).

Red-4 2.20 would be in the 4th drawer (Red-4), in the 2nd row, 20th [Slot](#) in the drawer Red-4.

Promedia Drawer Layout

Promedia Drawers come in several different sizes and drawer configurations. This page is based around drawers that have individual Slots dedicated to holding one Volume each, whether cased or not.

From:
<https://rtfm.tapetrack.com/> - **TapeTrack Documentation**

Permanent link:
https://rtfm.tapetrack.com/technote/promedia_labels?rev=1636496923

Last update: **2025/01/21 22:07**

