Labelling Promedia Drawers

Labelling the Promedia Drawers in the following manner allows the configuration of TapeTrack allowing operators to quickly locate a Sloted Volume, or locate a specified Slot to place a Volume in.

With TapeTrack Zones and Slotting information configured correctly, each Slot will be referenced with the Zone-ID, Drawer row and Slot number in the format of Zone-ID Row#.Slot#.

For example **Red-2 3.16** would be in the 2nd drawer (Red-2), in the third row, 16th Slot.

Red-4 2.20 would be in the 4th drawer (Red-4), in the 2nd row, 20th Slot in the drawer Red-4.

Promedia Drawer Layout

Promedia Drawers come in several different sizes and drawer configurations. This page is based around drawers that have individual Slots dedicated to holding one Volume each, whether cased or not.

From:

https://rtfm.tapetrack.com/ - TapeTrack Documentation

Permanent link:

https://rtfm.tapetrack.com/technote/promedia labels?rev=1636496923

Last update: 2025/01/21 22:07

