Scan To Slot Via Special Operations Slot Allocation

To use the Allocate to Slot via the Special Operations menu, the selected Volumes must be either moving to or confirmed in a Repository with allocated Zones.

While this process can be used any number of Volumes, you can only select a Slot for one Volume at a time so this process is best suited to one or low number of Volumes. To Slot multiple Volumes more efficiently see Scan to assign using Repository Usage Map.

Select Volumes To Slot

Select Slots

Refresh Display

From:

https://rtfm.tapetrack.com/ - TapeTrack Documentation

Permanent link:

https://rtfm.tapetrack.com/technote/scan_to_assign_slot_allocate?rev=1628644250

Last update: 2025/01/21 22:07

