

# Scan To Slot Via Special Operations Slot Allocation

To use the Allocate to Slot via the Special Operations menu, the selected Volumes must be either moving to or confirmed in a Repository with allocated Zones.

While this process can be used any number of Volumes, you can only select a Slot for one Volume at a time so this process is best suited to one or low number of Volumes. To Slot multiple Volumes more efficiently see [Scan to assign using Repository Usage Map](#).

## Select Volumes To Slot

## Select Slots

## Refresh Display

From:  
<https://rtfm.tapetrack.com/> - **TapeTrack Documentation**

Permanent link:  
[https://rtfm.tapetrack.com/technote/scan\\_to\\_assign\\_slot\\_allocate?rev=1628644250](https://rtfm.tapetrack.com/technote/scan_to_assign_slot_allocate?rev=1628644250)

Last update: **2025/01/21 22:07**

